



LEADERS IN INFORMATIONAL  
WATER & SOIL TEST KITS



1-844-TEST-H2O  
1-844-837-8420

## Well-Scan Sample Collection Instructions

Read All Instructions Before Collecting Samples

Do Not Collect Samples on Friday, Saturday, Sunday or Before a National Holiday



1. Remove gel pack from kit and freeze for at least eight hours before collecting your sample.
2. It is best to collect your sample as late in the day as possible and ship the sample to the lab the same day it is collected.
3. The sample needs to be received at the lab within 30 hours of sample collection, so you should ship your samples via overnight delivery. It's best to make these arrangements before

collecting your sample.

4. Determine the best site to collect your sample. Most often this would be your kitchen tap. Sometime you may be experiencing water quality issues in a particular area of the home, so you should collect at a tap closest to the problem.



5. Remove any aerator or filtering device from tap you are about to sample.
6. Turn water on and let run until noticeably colder about 5-10 minutes.
7. Fill large square plastic bottle to the neck and cap securely and return to the cooler.
8. Turn water off and disinfect the spigot using an alcohol swab.

9. Turn on water again and run a few moments until noticeably colder, this make take less time than before.



10. Remove protective seal from bacteria bottle, do not use if seal has been broken. Remove cap and hold while collecting the sample. Take care not to touch the inside of the cap/bottle or set it down, to avoid possible contamination. Fill the sample to the fill line, taking care not to overfill. Cap securely immediately after filling.
11. Return bottle to the cooler and add frozen gel pack.

12. Fill out the paperwork, making sure to include the date, time and sample ID. **The sample identification can be an address or a specific location like kitchen sink, bathroom tap.** The sample identification will be helpful when comparing your results from year to year.

13. Place paperwork in Ziploc bag and return to cooler.

14. Tape up box and ship to: TestMyWater.info  
812 Silvia Street  
Ewing, NJ 08628

Thank you for choosing TestMyWater.info  
1-884-837-8420  
info@testmywater.info